

# Quantumfield Guide



Laser Protocols

<b>Quantum Lasers</b>	<b>4</b>
<b>Quantum Laser – What is it?</b>	<b>4</b>
<b>How this works</b>	<b>4</b>
<b>Uses and Benefits</b>	<b>4</b>
<b>Unique Approach to Lasers and Healing</b>	<b>5</b>
<b>What Makes These Lasers Unique</b>	<b>5</b>
<b>The Quantum Laser Approach</b>	<b>6</b>
<b>Quantum Laser and Pulsar Laser Specifications</b>	<b>6</b>
<b>405 VioPulsar Probe</b>	<b>7</b>
<b>650 Red Laser Pulsar Probe</b>	<b>7</b>
<b>780 Infra-Red Laser Pulsar Probe</b>	<b>7</b>
<b>Quantum Laser Protocols</b>	<b>8</b>
<b>Introduction</b>	<b>8</b>
<b>The Unified Field Protocol</b>	<b>8</b>
<b>Cranial Sacral Wave</b>	<b>8</b>
<b>The Glands or Sacred Points</b>	<b>9</b>
<b>The Alchemical Nature</b>	<b>9</b>
<b>The Quantum Anatomy</b>	<b>9</b>
<b>Quantum Rejuvenation and Anti Aging</b>	<b>10</b>

<b>Pleomorphic Somatidian Life Cycle</b>	<b>11</b>
<b>Overview</b>	<b>11</b>
<b>Operational Procedures</b>	<b>17</b>
<b>Probe Operational Procedures</b>	<b>17</b>
<b>Session Length</b>	<b>17</b>
<b>Getting Started!</b>	<b>19</b>
<b>ESSENTIALS Section:</b>	<b>20</b>
<b>Beginning Laser Setting: Unwinding</b>	<b>21</b>
<b>Session:</b>	<b>21</b>
<b>Regenerative Uses</b>	<b>23</b>
<b>Regenerative Uses for the Pulsar Probes</b>	<b>24</b>
<b>Selecting A Setting</b>	<b>25</b>
<b>Settings</b>	<b>26</b>
<b>To Edit Settings:</b>	<b>27</b>
<b>Battery Instructions:</b>	<b>27</b>

# Quantum Lasers

## Quantum Laser – What is it?

The Quantum Laser is cold laser or low level laser system. This class of lasers is also often referred to as a soft laser or healing laser and refers to very low powered lasers which have been designed to gently rejuvenate cells. Low level lasers are totally safe and work on a completely regenerative basis.

## How this works

Cold lasers rejuvenate the cells by gently donating energy in the form of photons, which then increases the energy molecule of the cell, called ATP, which in turn stimulates the cell to grow more efficiently. This works in the same way that a garden works. Sunlight donates photons which in turn energizes the plants to grow bigger and healthier fruits and vegetables.

## Uses and Benefits

The Quantum Laser has been designed for the lay person and health practitioner alike. Anyone who wishes to use structured coherent light to rejuvenate cells will find great benefit from our line of cold lasers.

We specifically condition our lasers to not only rejuvenate the cell, but also to unwind and clear cell memory, to activate latent quantum biological processes. We do this with advanced patented scalarwave technology and unique violet crystal diodes.

Our systems are used by many health practitioners such as chiropractors, acupuncturists, naturopaths, biological dentists, physical therapist, sports teams and of course lots and lots of lay people around the planet.

They are used for injuries and pain relief, stress and dis ease, organ detox and balancing, cranial and glandular unwinding, acupuncture and energy healing, lymphatic wellness and rejuvenation and quantum activation and shift to the unified field.

The Quantum Laser and probes are easy to use and designed for everyone to be able to apply easily. We are not doctors and the Scalar Wave Protocols are not designed to treat or cure medical conditions. Rather this approach is based on the understanding that the body is fully capable of rejuvenating and finding solutions to its injuries and dis ease. Especially when its key quantum systems are relaxed and at ease, and when there is plenty of coherent energy available, either from lasers or the universe of which we are all connected.

## **Unique Approach to Lasers and Healing**

Basically it is our view that the universe is filled with coherent neutral energy, or space. When people become stressed and tense they hold onto cellular memory, which renders the body incoherent. When the body is incoherent it no longer taps into the vast reservoir of space and neutral energy.

Our approach to lasers is to unwind and clear this cellular memory so that the body can relax and unwind and open back up to the abundant energy of the universe. When it does this injuries, stress and dis ease disappear. As does energy when it flows and is abundant is vital and full of life.

So in summary our approach is not fix or heal the injuries or the “dis” “ease”, rather we seek to unwind the cell memory which tells those cells and body systems to remain closed off to the universe of energy. We unwind this stress and tension in much the same way that an old tape recorder is erased by simply reversing the polarity. The Quantum Lasers use a revolutionary scalarwave technology to achieve this ability to easily delete the cell memory.

## **What Makes These Lasers Unique**

All low level lasers donate energy in order to rejuvenate the cell. What makes our soft lasers truly unique is the ability to clear the cell memory. Cell memory is the stress or tension pattern underneath all injury and dis ease. And our proprietary technology activates scalar waves which neutralizes all polarity and shifts the cell memory. This quantum laser approach can be applied to any injuries, stress, tension, dis ease, or to activate the key body systems such as the glands, organs, chakras, meridians or any of the cells which make up the body.

## The Quantum Laser Approach

Basically the Quantum Laser approach enables the user to regenerate cells and to clear stress and tension pattern of cells in the body. It comes with built in protocols for all body systems where there is pain and a variety of other uses. The Quantum Laser Unwinding approach, quantum protocols, detoxifying qualities, regenerative protocols, and protocols for dealing with emotions, stress relief and the activation of consciousness, all make this laser a uniquely subtle experience with overall power and strength to see quick results. Some of the benefits Include:

- ◆ Pain Relief
- ◆ Circulation
- ◆ Arthritis
- ◆ Yogic Benefit
- ◆ Wellbeing
- ◆ Muscle Relaxation
- ◆ Stiffness Relief/Unwinding Effect
- ◆ Inflammation
- ◆ Stress and Tension
- ◆ Emotional Clearing

## Quantum Laser and Pulsar Laser Specifications

- ◆ 8 780 5mw infrared laser diodes
- ◆ 8 650 5mw red laser diodes
- ◆ 20 Violet 5mw spectrum LED diodes
- ◆ Inputs Digital Display Read Out
- ◆ Most Subtle Laser Technology
- ◆ Digital Processor
- ◆ Creates Rounded Sign Waves
- ◆ For Subtle Effect
- ◆ Phase Shift Capability Input within the Laser Beam Scalar Wave Variation

## 405 VioPulsar Probe

- ◆ 1 100 mw Violet Laser Diode
- ◆ 4 5mw Violet Spectrum Led Diodes
- ◆ Phase Shift Scalar Wave
- ◆ Unique Violet Spectrum Diode
- ◆ Diode holds more cellular Information than other colors
- ◆ More Subtle/Quantum Spectrum
- ◆ Shifts the Body into Unity State

## 650 Red Laser Pulsar Probe

- ◆ 1 100mw Red Laser Diode
- ◆ 4 5mw Violet Spectrum Led Diodes
- ◆ Full Phase Shift Scalar Wave
- ◆ Red Laser Diode for Soft to Medium Tissues

## 780 Infra-Red Laser Pulsar Probe

- ◆ 1 100mw Infra-Red Laser Diode
- ◆ 4 5mw Violet Spectrum Led Diodes
- ◆ Full Phase Shift Scalar Wave
- ◆ Infra-Red Laser Diode for Medium to Hard Tissue

# Quantum Laser Protocols

## Introduction

The Quantum Laser training process is based on a quantum understanding of the body called the Unified Field Protocol. This revolutionary low level laser training approach is based on the understanding that cells love to grow and are very good at doing this.

Anyone who tells you that the body does not regenerate is not telling you the truth, the only thing that the body does is to constantly regenerate cells. The Quantum Laser approach is based on using lasers to donate energy to regenerate the cell.

All of our protocols work on the understanding that in order to assist the cells and the body to do this work of regenerating cells we have to first unwind or clear the cellular energy. This cell energy is what limits new cells from the full state of energy and the quantum state wellness.

## The Unified Field Protocol

The Unified Field Protocol is a proprietary process for unwinding cells, organs, glands, meridians and chakras into the unified field. Basically all cells and body systems are designed to be relaxed and at ease in neutrality. The universe is 99 % neutral energy or space and when we relax and unwind the polarity, or cell memory, the cell opens up to this unified field of energy.

Our approach is to restore harmony to the field which is the underlying and essential nature of the body. In fact, the body is really a dimensional contraction of the field. When we understand the nature and operation of this dynamic relationship our work as stewards of health is to simply unwind, clear, detoxify, regenerate or activate the system.

## Cranial Sacral Wave

The protocols work with a systematic understanding of osteopathic principles in order to unwind the sympathetic nervous system and to reintegrate this into the parasympathetic nervous system which is where yogis and mystics reside.



When we are relaxed and at ease and the cranial sacral wave front is unwound, it shifts or neutralizes into the quantum state often called the still point. Yogis and Mystics refer to this place as the space between each breath. The Unified Field Protocol aims to return the cellular field or body to this quantum state of being neutral which naturally plugs us into the vast and neutral field of energy which makes up this universe.

## **The Glands or Sacred Points**

The glands or sacred points are key in the Unified Field Process. In fact the Unified Field Protocol centers around the glands and the understanding that the glands govern the cranial sacral wave front of the body and also condition everything from chakra plexuses to the nervous system and the cells of the body.

Our view is that the glands work through a principle called entrainment which simply means that they are conditioned by the people they have been around such as family. This is the same principle that happens when women live in a house together and their monthly cycles come together. The Unified Field Protocol unwinds this glandular entrainment in a systematic fashion which unwinds the cranial sacral system and all cells in the body simultaneously.

## **The Alchemical Nature**

As the adrenals are unwound the sympathetic nervous system unwinds and the alchemical and higher dimensional nature of the glands activates returning the body to a more coherent state of being. This process also activates key neurotransmitters and sacred substances such as oxytocin (released in breast feeding and orgasms) which further helps to free the body of stress and tension and induce more enjoyable states of being. The key is to clear out the excess adrenaline from the body and to activate the quantum anatomy and subtle hormones responsible for unity plane consciousness

## **The Quantum Anatomy**

As the body dissolves into neutrality, the subtle quantum anatomy activates. The basic protocol or sequence of the Unified Field Protocol is that the DNA induces the glands, which then conditions the cranial sacral wave condition, which in turn activates the plexuses of energy, the organs, and all the cells in a given area.

The Protocol also works great on all injuries, disease, stress and tensions by shifting the cell memory, or contraction of energy. This approach works great for everything from the smallest cells and DNA to the larger organs and even the emotions and entire chakra field of energy.

It does not matter what you are working on as the Quantum Laser works to shift or neutralize polarity into the neutrality. This phenomenon occurs in all phenomena from the particle/ wave, to this planet/star, to the body/unified field.

Remember our goal is to recognize that each of these are, in fact, one phenomena and conditioned when we approach them with a complete understanding of polarity, neutrality and scalar wave physics.

## **Quantum Rejuvenation and Anti Aging**

Adrenaline is the biggest epidemic on the Planet and is responsible for everything from chronic pain and stress, high blood pressure, diabetes, disease and most of all aging.

Studies have shown that when the adrenals are constantly releasing adrenaline and cortisol this shuts down the other glands including the thymus or heart chakra gland responsible for immunity and anti aging. What happens is that whenever adrenaline is released the cells get the signal to turn on then off. As this happens they shift polarity and subsequently lose their charge, their oxygen and energy and aging occurs.

This process shortens the telomeres or scalar antenna on the DNA which as they get shorter limits the ability to grow new cells, until life cannot go on in this form. Fortunately scientists have recently discovered an enzyme called the Telomerase enzyme which increases the length of the telomere thereby increasing the anti aging potential and giving hope for the potential of longer life cycles on this planet.

It is our view that both Scalar Waves and the violet laser spectrum optimizes and enhances this quantum rejuvenation effect of activating Telomerase and the Telomere length. Our protocols also aim to activate stem cells, key neurotransmitters and sacred compounds of the body, and further refine and activate the genetic sequencing capabilities of the quantum field.

In addition to all of the amazing health and wellness applications of the Quantum Laser our favorite applications have to do with the subtle energy and shift toward unity plane consciousness. Many who use our lasers start off by using lasers for pain and injuries, and then

expand their uses to the subtler work of dealing with emotions, the chakras, quantum healing, and unity plane consciousness.

## **Pleomorphic Somatidian Life Cycle**

The Quantum Laser Protocols are based on a regenerative principle of wellness. The pleomorphic model views blood as living (not sterile like the western model) and the goal as returning the blood to a non pathogenetic state. Our goal is not to destroy the germ, rather we seek to unwind the terrain and return the somatid to a non toxic state.

The Somatid is the smallest unit of life energy. These are the little points of light which can be seen when viewing blood under a darkfield microscope. These points of light have been described as indestructible and eternal.

The Somatidian life cycle has been mapped out using very high powered microscopes and can be viewed, as in either regenerative states or degenerative states. In the latter the somatids actually cycle through different organistic phases, which become the germs which the germ theory tries to destroy.

Once again our aim is to return these points of light to a non toxic regenerative state of energy. When in this state there is optimal health and an eternal state of well being.


## **Overview**

The Quantum Laser comes with preset phase shift scalar waves for dealing with all the cells of the body and all of the systems of the body.


The Unified Field Protocol is not based on germ theory or the attempt to fix or cure. It is based on the principle and understanding of pleomorphic principle that the body finds solutions to situations when it is relaxed and at ease, and when all of it's systems are up and running, and working together. The body is all light and the Unified Field Protocol is designed to return all cells, organs and the entire body/field to coherence, natural wave structure and return to the Unity Plane consciousness.




Unified Field




## Unwinding Protocol



Organs



# Chakras



Nadis



# Operational Procedures

<b>STAR</b>	Press Once to Turn On. Hold Down to Turn Off.
<b>Arrow Right</b>	To Select an Item Moves Curser to the Right
<b>Arrow Left</b>	Goes Back to Previous Menu Saves Design or Edit Settings Moves Curser to the Left
<b>Arrow Down</b>	Moves Curser Down
<b>Arrow Up</b>	Moves Curser Up

## Probe Operational Procedures

Warning: Do not look directly into Pulsar Probe lasers. Lasers not intended for use directly on or near the Eyes.

Plug the Probe wire into center plug outlet of the Quantum Laser. Now operate choosing whichever phase shift setting you want. With the probe wire plugged into the main unit, the device will automatically filter the power and light through the probe rather than the main unit. The probes also run on 4 minute cycles. When the probes power off, press the right arrow several times to operate again on the same setting, or press the left arrow button to go back and chose a new setting. Or press the star button to turn off.

## Session Length

The Quantum Laser runs for 1 cycle, which is 4 minutes, then turns off. To restart the session, arrow over 2 to 3 times and it will start again for another cycle. When working on an injury, you may do 2 to 4 cycles which translates into 8 to 16 minutes. That is a good length for a session. After a rest you may go longer doing another 1 to 2 cycles. You can do this several times a day

if you are really in pain or trying to turn something around quickly. Usually after 10 to 20 minutes, you should feel a difference. If you do not feel a difference after several minutes on one cycle, try a different setting. Some situations need to unwind or relax at first, while others may need a pain relief setting. You may use your intuition when picking a setting. Most are very obvious as they work on specific areas of the body and field.

# Getting Started!

To start your laser simply press the STAR button.

A welcome screen appears and then shifts to a menu screen which says:

LAUNCH DESIGN

SETTINGS WELCOME

By scrolling with the arrow buttons you can navigate up or down. Pressing the right arrow will open up a sub menu and/ or start a program/function, while pressing the left arrow will return back to the previous menu and/or turn off a given program/function.

The Launch Section is the place where you can choose specific preset programs.

The Design Section gives you the ability to design, try and save your own frequency programs.

The Settings Section gives you the ability to change the setting on your laser.

The Welcome Button simply refreshes the Welcome Screen and then returns you to this same menu screen.

If you scroll right to the Launch Section you will come to another menu screen which says:

ESSENTIALS

# ESSENTIALS Section:

**Unwinding** unwinds cell memory, stress and tension from the body, cells and energy. It is great to use on old injuries and on any situation which has stress, tension and contraction. The basic gist of this program is the dissolving of contracted cellular memory and the return to the natural state of ease.

**Quantum** is about shifting energy back into the quantum state. It works in a similar fashion to the Unwinding Program only it does so on the energetic level. It is a very popular program which works great on all cells, organs, chakras and anywhere that is tense and contracted.

**Relieve** is a favorite for pain and injuries. It works to dissolve and relieve stress, tension, injuries and pain.

**Pleo** is a detox and immune activation program. It is based on research of the pleomorphic somatidian life cycle – and in particular the regenerative state of the somatid. Research has shown that these tiny point of light within the blood go through a number of different life cycles some indicative of dis ease while others reflective of health and wellness. Our goal here is to unwind, normalize and activate the somatid to its regenerative state.

**Activate** is based on higher dimensional quantum DNA physics and the principle of subtle regenerative coherence. This program is great for injuries, pain, and any cellular regeneration.

**Rejuv** is designed to activate the lymph system. In particular it is to induce the electrons into a high spin state. When electron shift into a high spin state they phase shift into an ether state such as the lymph, cranial sacral, brain, sexual and glandular substances. This program is really a higher dimensional wellness program to induce this high spin alchemical state.

**Soothe** is a very soothing program for injuries, pain and stress. It works great by dissolving tension.

**Unity** is a beautiful program for neutralizing anything and everything - and shifting into unity. It works wonderfully through the violet laser light in processing in-form-ation into the Quantum State.

**Vital** is about activating the inner cellular energy. It is great for activating the immune response, activating oxygen and electrons, detoxification and clarifying energy, regeneration and wellness.

**Cohere** Program brings structure to the body and field. It is about organizing, processing and actualizing energy in-form-ation. It normalizes all cells, body systems and chakras into the one integrated field space. The body is designed in and of coherence and all health and or dis ease is simply the reflection of how coherent we are and how efficiently we can process information (in-form-ation).

**Wave** is for dissolving all cells and phenomena into one wave essence and then allowing that wave to expand into consciousness.

The CUSTOM section enables the user to be able to try, program and save their own protocols.

## **Beginning Laser Setting: Unwinding**

Ideal Position: Face Down on Massage Table, or Bed, etc... Unwinding Session Steps:

1. Place the Laser on the Unwinding/Quantum Setting over one lower back around the adrenals (we access them right above each kidney). Does not matter which side of the lower back you start the session with. Run for 4 minutes per side.
2. Place the Laser on the Unwinding/Quantum Setting on the sacrum (balanced on the tail bone). Run for 4 minutes on the Sacrum.
3. Place the Laser on the Unwinding/Quantum Setting on the back of the neck (we often turn people over on the table and rest the laser on it's side nestled right into the side of the neck.). Run for 4 minutes on the Back of the Neck.
4. Turn person over and Place the Laser on the Quantum Setting on the Thymus Gland (located on the lower/center of chest). Run for 4 minutes.

## **Session:**

Some people can unwind for quite a while. We usually do at least 1 or 2 minutes on each adrenal and can do many more if the person is unwinding nicely. To see if someone is

unwinding nicely, we set the laser on the first adrenal and watch for deeper breaths and fuller breaths as the person is lying face down on the table.

If someone has a lot of emotional stuff going on and is not immediately opening up to nice big full breaths, then we switch the setting to the Quantum Setting instead of the Unwinding Setting in the beginning which helps to clear out energetic holding patterns.

Then we move the laser to the second adrenal and do the same setting usually on that adrenal. After a couple of minutes on the Quantum Setting, we can switch back to the Unwinding Setting to work on each adrenal for a minute or two for the nervous system as well. Usually after a quick clearing of the energetic field, the nervous system is more flexible and ready to unwind too.

Next we move the laser down to the sacrum to unwind any holding patterns and release the sacrum. Balancing the laser on the tail bone, we start with the Unwinding setting and can do the Quantum Setting a bit too. We Run the laser on the sacrum for 2 to 4 minutes usually.

Next we put the laser near the back of the cranium on the where the head meets the spine. This unwinds the occipital ridge area which relaxes and opens the vagus nerve. This is a very important nerve and area as it also opens up the cranium, unwinds the hypothalamus (which governs the endocrine system) and relaxes the TMJ.

Lastly, we set the laser on the thymus gland (right on the center of the chest or breast bone) on the Quantum to strengthen the immune system, open up the heart chakra and lift emotions. It is a nice way to finalize a session.

This Unwinding Protocol is the underlying essence of much of the Quantum Laser approach. It is our view that the nervous system and cranial system holds stress and tension in the form of cellular memory. In order to systematically unwind and clear this stress and tension this protocol goes straight to the cranial sacral system, the glands and the subtle unified field.

This is very important if people have had a lot of stress in life and in particular if someone did not breast feed when they were an infant. Many folks who struggle with chronic pain, immune disorders such as chronic fatigue, ongoing stress, migraines, emotional difficulties, chemical sensitivities, etc. never Breast feed.

Breast feeding pulls the TMJ forward during sucking which relaxes the occipital ridge, which then releases the vagus nerve. The vagus nerve is the key parasympathetic nerve in the body

and pretty much does everything in your body like breathing, digestion, circulation, organ functions, etc.

For people who don't breast feed they carry this stress instead of unwinding it as they move through life. A little child who bumps into something will then cry, and then breast feed. This process relaxes the vagus nerve and the mothers milk also has a lot of oxytocin which signals the hypothalamus to come on and relax and dissolve the stress.

After doing this for a year or two the child is well versed in how to unwind stress, from a neuro endocrine level and they move through life handling stress. For people who never breast feed they tend to hold their stress and then go on to the next thing, until they reach their middle ages and suddenly have a major health challenge like chronic fatigue.

So basically we all could use some unwinding and clearing of stress, and most importantly, to reset our nervous system and glands so that we can be relaxed and at ease. As this is how the body and cells do their best at regenerating wellness.

## Regenerative Uses

### Injuries:

---

**For Pain Relief/Healing:**

Essentials: Relieve, Soothe, Activate

Location: Directly over Injury or Pain

---

**Release Injury From Trauma:**

Essentials: Unwinding, Quantum

Location: Directly over Injury or Trauma

### Immune System/Infections

**To Boost Immune Setting:**

Essentials: Vital, Unwinding, Quantum

Location: Thymus Gland, Plexus, Chakra,  
Meridian Area

---

---

**To Clear Micro Biological Infections/  
Arthritis**

Essentials: Pleo, Vital

Location: over affected area

## **Regenerative Uses for the Pulsar Probes**

### **Red 650 Pulsar Probe**

The 650 Pulsar is great for surface level cells and tissue such as soft tissue, muscle, nerve, meridians, skin, gums, etc. It is also great for stimulating blood, activating key enzymes like the cytochromase oxidase enzymes responsible for mitochondrial ATP production, increasing oxygen and electrons, pushing toxins into the lymph and for acupuncture meridian protocols. This is a great all purpose super useful wavelength and soft healing cold laser probe.

### **780 Infra-Red Probe**

The 780 Pulsar is designed for deeper seated cells such as bone, tendon, cartilage, disc, teeth, etc. The length of the wave sends the photons right through the surface cells and into the deeper more intrinsic cells. This wavelength is very popular for bone and teeth regenerative work. It is also used extensively in the stop smoking and weight loss protocols. We infuse our infra red Pulsar Lasers with violet spectrum light and Scalar Wave technology to help make them very subtle and gentle.

### **405 Blue-Violet Pulsar Probe**

Our soft laser approach is to use violet laser diode technology to activate the subtle quantum anatomy. In particular for regenerative purposes such as soft and hard tissue injuries, DNA activation, anti aging, facial rejuvenation, pleomorphic detoxication, subtle energy alchemy, and a vast variety of other low level laser therapy applications. Our view is that the violet spectrum induces the bodies quantum regenerative capabilities to kick on. Studies show that violet lasers are being use for activating the Telomerase enzyme and the Telomeres, activating stem cells and DNA. Our view is the violet laser opens the higher glands and activates subtle body anatomy.



## Selecting A Setting

The protocols that we offer are really quantum in nature and design, and center around either the body and the subtle anatomy or unified field. As we share our reference or experience, we encourage you the user to try things and play around with different possibilities. Allow your experience to guide you.

Our approach to research and protocol design, is to try lots of things and to see which ones feel good and also yield good results. Our main goal in all our protocols is to unwind stress and tension and to allow the body to relax. So for us it is essential for our programs to feel relaxing and shift or neutralize cell memory.

Our approach to using lasers or any other wellness protocols is to try something and notice the results. We often check in with our clients after a cycle or two of using the laser to see how it feels. If it is feeling good and like it is helping we continue with that program. If not we will then try a different program setting To give you an example, someone comes in with a pain in their knee. Our first approach is to unwind and clear gcell memory with either the Unwinding Program or the Quantum Program. Often these alone will relieve the pain and help to make the area feel better and also rejuvenate cells. In some cases we might then go to the Relieve Program which is great for pain relief and also balances the neurotransmitter sequencing with the brain to change the pain signaling into something different like, say, natural opiate (pain relief) signaling. Or the Soothe Program which stimulates serotonin is also great for pain. If we don't see the client unwinding or relaxing, or taking nice breaths and smiling then we might check in with them and see how the protocol is feeling. We are looking for them to feel some improvement and to feel good. If we don't get any such indication, we often try another program.

It is a little like watching the TV. If you don't like a channel that you are on what do you do. You try another one. Well this works the same way. As you use the Quantum Laser you will develop a sense of when it is having this unwinding effect. Eventually you will even be able to tell when it the right program on those who you work on with the laser.

For us there is a definitive response that a persons nervous system has when it is dropping cell memory and relaxing and that is what we are seeking.

For some who come in with the very same knee pain, we say well this mode isn't yielding response that we are looking for. And then we might try another program like the Pleo Program and then suddenly we see them breathe more fully and smile and so we keep going with that setting. The key and goal here is to unwind. When this happens it is clear and we would then view this issue as potentially having to do with a low grade infection, toxicity or something other than the an injury or physical stress alone.

Or we might try the Activate Program, and if that yielded good results we might consider that there was poor circulation or lack of blood. This could limit the flow of oxygen and electrons which could diminish the rebuilding of health bone in a knee.

We are not doctors and we don't share these as diagnosis, rather we are sharing how we approach laser rejuvenation. Our goal is not to diagnose or to fix. We simply try different things until we find the program which has the result of enabling the persons own system to unwind and rejuvenate their own cells. Which the body is perfectly capable of doing and has been doing since it arrived here on this Planet.

Remember this is quantum and results can often be felt immediately. So check in and be flexible and try different programs. Most conditions will respond nicely to any of the settings and programs, as we condition them all with unique Scalar Waves and our Unified Field research.

## Settings

Battery Remaining:	1 bar is low, 4 bars if full
Session	Suggested at 4 min/ Default at 4 min (# of minutes for each cycle) 1-20
Sound	(0-1) 0=off 1=on
Contrast	Suggested at 11% Varying capability of contract on the display screen

---

Backlight	Suggested at 100% Varying capability of contrast on the display screen
Wave	(0-359) Scalar Wave degree shift

---

## To Edit Settings:

1. Arrow to Right to Edit
2. To Edit Mode Type in Requested Number
3. Left Arrow Back Once to Save

## Battery Instructions:

The Unit may come with some charge in the battery already. You DO NOT need to charge the unit before using it for the first time. For optimal battery performance, unplug after it completes charge.

## ***QiFi Purifier Guidelines***

The right side knob is for selecting the setting. There are 6 settings. The left side knob is for fine tuning the setting.

There are 10 clicks to go full circle with the right knob.

The settings (first 6 clicks) go from 9:00 to 3:00 of a clock (there is no 12 noon position). 9:00 is the 1st of the 6 programmed settings.

The last 4 clicks that do not have a frequency assigned to them are just pure clean quantum field generated light. They still can be used since they have the crystal oscillators that help create the quantum field to harmonize the EMFs.

### Settings & Descriptions

**Unwinding** Unwinds cell memory, stress and tension from the body, cells and energy. It is great to use on old injuries and on any situation which has stress, tension and contraction. The basic gist of this program is the dissolving of contracted cellular memory and the return to the natural state of ease.

**Quantum** is about shifting energy back into the quantum state. It works in a similar fashion to the Unwinding Program only it does so on the energetic level. It is a very popular program which works great on all cells, organs, chakras and anywhere that is tense and contracted.

**Activate** is based on higher dimensional quantum DNA physics and the principle of subtle regenerative coherence. This program is great for injuries, pain, and any cellular regeneration.

**Unity** is a beautiful program for neutralizing anything and everything - and shifting into unity. It works wonderfully through the violet laser light in processing in-form-ation into the Quantum State.

**Cohere Program** brings structure to the body and field. It is about organizing, processing and actualizing energy in-form-ation. It normalizes all cells, body systems and chakras into the one integrated field space. The body is designed in and of coherence and all health and or dis ease is simply the reflection of how coherent we are and how efficiently we can process information (in-form-ation).

**Wave** is for dissolving all cells and phenomena into one wave essence and then allowing that wave to expand into consciousness.

The training around using the Qifi Purifier is to select a setting and just see how you feel and then keep selecting until you land on one that gives you the feeling you desire in this present moment. You can run it for any length of time, even 24/7. The range for the Qifi Purifier is more than 3/4 of a mile.

The fine tuning offers more of the energy (turning the dial up or to the right) or less of the energy (turning the dial down or to the left). This dial has a stop whereas the right dial clicks all the way around. We (the company staff) tend to keep the knob in the center.

## ***All-Violet Lotus Laser for Subtle Energy Balancing & Extraordinary Healing***

The Lotus Laser operating settings are exactly the same as the main ScalarWave Laser device.

Similar in design to the main ScalarWave Laser unit, the Lotus Laser is the first of its kind all-violet, digital, hand-held, portable laser for anti-aging, rejuvenation, skin care, and pain relief.

Violet laser light is highly compatible with subtle energy healing and vibrational medicine methods, such as sound therapy, color therapy, chakra balancing, acupuncture and acupressure, and many others.

The Lotus Laser is a landmark laser unlike any that has ever come before it—a quantum therapeutic healing device. It features 16 true violet laser diodes and 20 super violet LEDs, all of which emit therapeutic violet wavelength energy in the violet spectrum.

R&D Research shows that the violet wavelength has special properties that are being used to:

- Supports DNA
- Increases length of telomeres
- Supports stem cells
- Balances subtle energy, positively affecting everything from the emotions to enzyme functions to the energy fields surrounding the body, and more

